

# Contents

An introduction to the custom preprocessing pipeline built for surfaced-based analysis

1. Motivation – Why surface?
2. Overall scheme – Covering multiple softwares (e.g., FSL, AFNI, ANTs, Freesurfer)
3. Pipeline details – Comparison with existing preprocessing pipelines
4. Practical applications – Not limited to surface analysis



### 3. Pipeline details – Comparison with existing preprocessing pipelines

#### Overall scheme of this pipeline

A) Make directories for preprocessed data (r1)

B) Convert DICOM images to NIFTI format (using `dicm2nii.m`, which was adapted from <https://github.com/xiangruili/dicm2nii>) (r2-r4)

C) Basic environmental setup for preprocessing (s1)

D) Preprocessing structural data

For surface-based preprocessing:

1. Use `recon-all` (Freesurfer) to correct bias-field, extract brain tissue, reconstruct structural images to cortical surface, and do anatomical segmentation (s2)
2. Use `ciftify_recon_all` (CIFTIFY, Dickie et al., 2019) to normalize structural images to MNI space, and resample native surface images onto Conte69 164k and 32k CIFTI surface (Van Essen et al., 2012) using MNI normalization parameters (s3)


For Volume-based preprocessing:

1. Use `antsBrainExtraction` (ANTs) to correct bias-field and extract brain tissue (s4)
2. Use `antsRegistrationSyN` (ANTs) to normalize structural images to MNI space (s4)
3. Use `FAST` (FSL) to do anatomical segmentation (s4)

- Similar with the SPM-based COCOAN pipeline.
- For surface-based preprocessing, `recon-all` and `ciftify` was used to project structure data from volumetric space to CIFTI space.
- For volume-based preprocessing, ANTs were used.




### 3. Pipeline details – Comparison with existing preprocessing pipelines



<https://www.youtube.com/watch?v=zeFPx0fMXRQ>

- |                                       |  |
|---------------------------------------|--|
| 1. Motion Correction                  | 17. Tessellation   |
| 2. NU Intensity Correction            | 18. Orig Surface Smoothing                               |
| 3. Talairach                          | 19. Inflation  |
| 4. Normalization                      | 20. QSphere  |
| 5. Skull Strip                        | 21. Automatic Topology Fixer                             |
| 6. Automatic Subcortical Segmentation | 22. Final Surfaces                                       |
| 7. EM (GCA) Registration              | 23. Cortical Ribbon Mask                                 |
| 8. CA Normalize                       | 24. Spherical Inflation                                  |
| 9. CA Register                        | 25. Ipsilateral Surface Registration (Spherical Morph)   |
| 10. Remove neck                       | 26. Contralateral Surface Registration (Spherical Morph) |
| 11. EM Registration, with Skull       | 27. Average Curvature                                    |
| 12. CA Label                          | 28. Cortical Parcellation                                |
| 13. ASeg Stats                        | 29. Parcellation Statistics                              |
| 14. Normalization2                    |  |
| 15. WM Segmentation                   |  |
| 16. Cut/Fill                          |  |



[https://surfer.nmr.mgh.harvard.edu/fswiki/FsTutorial/OutputData\\_freeview](https://surfer.nmr.mgh.harvard.edu/fswiki/FsTutorial/OutputData_freeview)




### 3. Pipeline details – Comparison with existing preprocessing pipelines

#### cifti-for-your recon\_all outputs

##### ciftify\_recon\_all

- T1w/Native: The freesurfer "native" output meshes
- MNINonLinear/Native: The T1w/Native mesh warped to MNINonLinear
- MNINonLinear/fsaverage\_LR32k
  - the surface registered space used for fMRI and multi-modal analysis
  - This 32k mesh has approx 2mm vertex spacing
- MNINonLinear\_164k\_fs\_LR (in the MNINonLinear folder):
  - the surface registered space used for HCP's anatomical analysis
  - This 164k mesh has approx 0.9mm vertex spacing




| surface       | description  | useful for  |
|---------------|--|---|
| white         | The border between the gray matter and the white matter                      | defines the inside of the cortical ribbon                         |
| pial          | The border between the gray matter and the outside of the brain (pia matter) | defines the outside of the cortical ribbon                        |
| midthickness  | The midpoint between the white surface and pial surface                      | measuring distance and surface area                               |
| sphere        | The surface vertices as a sphere   | registration and resampling                                       |
| inflated      | The midthickness is blown up like a balloon                                  | visulation  |
| very-inflated | The midthickness is blown up like a balloon <i>more</i>                      | visualizations  |
| flat          | the cortex ripped apart  | visualizations (one less intuitive view of the entire hemisphere) |

[https://edickie.github.io/ciftify/#/03a\\_cifti-for-your\\_recon\\_all](https://edickie.github.io/ciftify/#/03a_cifti-for-your_recon_all)




### 3. Pipeline details – Comparison with existing preprocessing pipelines



### 3. Pipeline details – Comparison with existing preprocessing pipelines

SPM (cocoanlab) vs. ANTs



### 3. Pipeline details – Comparison with existing preprocessing pipelines

#### E) Preprocessing functional data

1. Use 3dTshift (AFNI) to do slice-timing correction (s5, if needed)
2. Use 3dvolreg (AFNI) to do motion-correction (s6)
3. Use topup/applytopup (FSL) to do distortion-correction (s7)
4. Use flirt (FSL) with BBR cost function to co-register functional images to structural images (s8)
5. For surface-based preprocessing, additionally use bbregister (Freesurfer) to refine co-registration (s8)
6. Use ICA-AROMA (Pruim et al., 2015) to remove motion-related signals (s9)
7. Use 3dTproject (AFNI) to remove nuisance signals (s10)
8. Use applywarp (FSL, for surface-based preprocessing) or antsApplyTransforms (ANTs, for volume-based preprocessing) to normalize functional images to MNI space
9. For volume-based preprocessing, use susan (FSL) to spatially smooth functional images
10. For surface-based preprocessing, use ciftify\_subject\_fmri (CIFTIFY, Dickie et al., 2019) to transform functional images to Conte69 32k CIFTI surface, and spatially smooth (based on both surface-based and volume-based smoothing) functional images

Note: 8 can precede 6 and 7. depending on your purpose.

- Similar with the SPM-based COCOAN pipeline.
- For surface-based preprocessing, fMRI volumes were projected to CIFTI surface!
- For volume-based preprocessing, outputs are almost same.







### 3. Pipeline details – Comparison with existing preprocessing pipelines

#### Motion correction: AFNI 3dvolreg


Realignment accuracy



Interpolation accuracy



Speed




Oakes et al, 2005





### 3. Pipeline details – Comparison with existing preprocessing pipelines

#### Coregistration: FSL FLIRT + BBR (Boundary-Based Registration)



Greeve and Fischl, 2009



### 3. Pipeline details – Comparison with existing preprocessing pipelines

#### Coregistration: FSL FLIRT + BBR (Boundary-Based Registration)

Pre-alignment by FSL FLIRT (cost function: correlation ratio)



Alignment by FSL FLIRT (cost function: bbr)



Refinement by freesurfer bbregister (cost function: bbr)


(if freesurfer reconstruction  
was done before)



### 3. Pipeline details – Comparison with existing preprocessing pipelines

#### Denoising: AFNI 3dTproject


‘one-step’ strategy is the best!



Lindquist et al., 2019

steps within an omnibus framework. For example, it is relatively straightforward to formulate a single linear model that simultaneously performs motion regression, nuisance regression, and temporal filtering. This is an approach advocated by Caballero-Gaudes and Reynolds (2017), and implemented in Analysis of Functional NeuroImages (AFNI) (Cox, 1996) **3dTproject**. In general, we believe that the

Instead of bpf -> Nuisance regression, combine them onto one matrix.



### 3. Pipeline details – Comparison with existing preprocessing pipelines

#### **Normalization:**

**FSL applywarp (surface-based)**

**or ANTs antsApplytransform (volume-based)**

Surface-based: Use pre-computed transformation parameters (FSL FLIRT+FNIRT) from ciftify\_recon\_all step.

Volume-based: Use pre-computed transformation parameters (ANTs antsRegistrationSyN) from normalization step.



### 3. Pipeline details – Comparison with existing preprocessing pipelines

#### **Smoothing: FSL SUSAN**

For volume-based analysis; Anatomically-constrained volumetric smoothing technique!



### 3. Pipeline details – Comparison with existing preprocessing pipelines

#### **Projection to surface and Smoothing: CIFTIFY**

For surface-based analysis; project fMRI volumes (in MNI space) to CIFTI surfaces (in Conte69 32k space) and perform 2D (cortex) + 3D (subcortex) smoothing!

